
8/96-Se-

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The device makes it possible to investigate the expansion of water in the temperature range from 0 °C to approx. 15 °C, whereby the height h of the water column in the capillary tube is a measure of the change in volume. By plotting a curve of $h = f(\vartheta)$, the anomaly of water can be demonstrated and a density maximum can be detected.

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1 Safety notes

When volume increases abruptly (indication of initial freezing of the water), **immediately** remove device from freezing mixture.

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2 Description, technical specifications

- ① Glass vessel
filling volume approx. 310 cm³
- ② Filler tube with tap (2.1)
diameter of hose bead: 8 mm
- ③ Capillary tube, 35 cm long, with mm scale
inner diameter: approx. 1.7 mm
- ④ Threaded tube for screw-type fasteners ⑤
- ⑤ Screw-type fasteners for threaded tube ④; with fittings (8 mm dia. and 1.5 mm dia.) for thermometer and temperature sensor respectively
- ⑥ Stirring rods, approx. 1.5 cm long

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Magnetic stirrer	666 845
Tubing, 7 mm inner dia.	667 180
Funnel	665 009
Stand base, V-shape, 20 cm	300 02
Stand rod, 44 cm	300 42
Clamp with ring	301 10
Coloring agents	666 609/666 572
Distilled water	
Dish soap	
Refrigerator	

or freezing mixture (approx. 450 g crushed ice, mixed with approx. 40 g table salt) in plastic trough (306 90).


Fig. 3

3.2 Experiment preparation

Insert stirring rods ⑥ into device via threaded tube ④; add some coloring powder and 1-2 drops of dish soap. Slide screw-type fastener ⑤ over thermometer or temperature sensor. Insert thermometer into the device through the threaded tube such that temperatures as low as 0 °C are clearly visible. Carefully tighten the screw-type fastener.

Connect funnel to the filler tube ② according to Fig. 2, add distilled water and open tap (2.1). Add water continually so that the funnel is full at all times while the device is being filled;

Remove air bubbles by swirling the device; when the water has reached the upper edge of the capillary, close tap (2.1) and remove tubing (if necessary, first pour remaining water out of funnel); by briefly opening the tap (2.1), allow the water level in the capillary to drop to a height of approx. 18 cm (for "refrigeration method" see below) or 30 cm (for "freezing mixture method").


Fig.4

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Recommended ways of obtaining necessary temperature range (cf. Fig. 4):

"Refrigeration method": The device, which has been filled with water, is allowed to cool overnight in the bottle compartment of the refrigerator; set thermostat such that the temperature here is approx. 1 °C (check beforehand).

Remove device from refrigerator immediately before experiment is commenced and place on magnetic stirrer.

"Freezing mixture method": Place device in freezing mixture prepared from 450 g chopped ice and 40 g table salt (mix thoroughly) in the plastic trough or similar vessel (Fig. 3).

Immediately afterwards, place device on magnetic stirrer.

3.3 Experiment procedure:

Turn on magnetic stirrer immediately and set to moderate speed (approx. 350 rpm).

Determine the water level h in the capillary as a function of the temperature (Fig. 4);

When freezing sets in (indicated by sudden increase in volume), immediately remove device from freezing mixture.