
7/97-Sf-

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Fig. 1

The Geiger counter (546 281) is used in introductory experiments on radioactivity (detection of α and β radiation; explaining the function of a pulse ionization counter tube).

The output pulses of the Geiger counter are either

- displayed using a Wulf electroscope (546 00) or
- made visible using an oscilloscope via the adapter for Geiger counter or made audible using a loudspeaker (e.g. 587 08) via a power amplifier (e.g. 522 61).

Literature:

Central Card-file: Radioactivity (598 691), in German

1 Safety notes

- Be careful when experimenting with high voltage! Do not touch voltage-bearing metal parts! Ground the housing of the Geiger counter!
- When using a radioactive preparation, observe all safety regulations which apply in your country and/or local jurisdiction.

2 Description, technical data

2.1 Geiger counter (546 281), Fig. 1

- ① Needle electrode, mounted in insulator (1.1), connected to 4-mm socket (1.2) for applying the positive pole of the high voltage source
Operating voltage: approx. 2.5 V - 3 kV
Pulse ionization range: < 100 V
- ② Metal housing (cathode) on 4-mm pin (2.1) for attaching to the Wulf electroscope (546 00) or the adapter for Geiger counter (546 38); with radiation aperture (2.2) (dia. 0.8 cm)
housing dimensions: 5 cm x 2.5 cm dia.
- ③ 10 plastic foils, vapor-deposited aluminum coating, approx. 0.005 mm thick, for covering the radiation aperture (dust protection)

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2.2 Adapter for Geiger counter (546 38)

- ④ Highly insulated 4-mm socket for plug pin (2.1) of Geiger counter (pulse input)
- ⑤ Highly insulated 4-mm pulse output socket for connecting the signal input of an oscilloscope or an amplifier
- ⑥ Ground socket

Housing dimensions: 13.5 cm x 7.5 cm x 5.5 cm

Fig. 2

3 Operation

Additionally required:

High voltage power supply 10 kV	521 70
Pulse indicator, e.g.	
Wulf electroscope	546 00
or	
Oscilloscope	e.g. 575 211
or	
AC/DC amplifier	522 61
with loudspeaker	e.g. 587 08

Recommended preparations

Radium preparation 3.3 kBq	559 59
Ra-226 preparation 3.3 kBq	559 430

Stand material

Carefully clean needle electrode ① of dust or similar using a damp leather rag or a flame (do not heat red-hot!). After inserting the clean needle, cover opening (2.2) with a foil (from ③).

If necessary, clean the insulators of sockets (1.2) as well as ④ and ⑤ using a lint-free cloth moistened with a little dish soap; dry the sockets with hot air.

Set up the apparatus as shown in Fig. 3 resp. Fig. 4.1 or Fig. 4.2. Connect the positive pole of the high-voltage source to the needle electrode ① via socket (1.2).

Setting up the pulse indicators:

Wulf electroscope: slide the cathode close to the strip.

Oscilloscope: vert. input: ac
vert. ampl.: 0.2 v cm^{-1} to 2 v cm^{-1}
time base: 1 ms cm^{-1} to 5 ms cm^{-1}

AC/DC amplifier: AC; gain x 1000

Setting the optimum high voltage:

Before placing the preparation in front of the radiation aperture, increase the high voltage very slowly to exactly the value at which spontaneous discharge occurs (indicated by constant deflection of the leaf toward the cathode of the Wulf electroscope as well as by numerous optical or acoustic pulses on the oscilloscope or speaker respectively). Then reduce the voltage by about 30 V - 50 V so that discharge can only occur when the air between the needle electrode and the housing of the Geiger counter is ionized by radiation.

Recommended distances between preparation and radiation aperture: approx. 1 cm to 3 cm for α and β radiation, 3 cm to 6 cm for β radiation;

To demonstrate β particles, hold a sheet of paper between preparation and radiation aperture to shield the α radiation.


Fig. 3

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Demonstrating α and β particles with the Wulf electroscopes (546 00)


Fig. 4.1

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Optically demonstrating α and β particles with an oscilloscope

Operating voltage in plateau of counter-tube characteristic: display of α and β particles by pulses of equal amplitude.

Operating voltage in proportional range: energy-analytical display (low pulses due β radiation, high pulses due to higher-energy α radiation).

Recommended for setting the high voltage in an extremely narrow proportional range: fine adjustment via 10-turn potentiometer 1 k Ω (577 93) at control input (a)/(b) (see inset drawing and Instruction Sheet 521 70).


Fig. 4.2

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Acoustically demonstrating α and β particles via AC/DC amplifier (522 61) and loudspeaker (587 08)